

Product Data Sheet

Product Name: PROTAC CDK6 Degradar 1

Cat. No.: GC80885

Chemical Properties

Cas. No.


Formula C₇₀H₈₉N₁₅O₈ M.Wt 1268.55

Solubility DMSO: 50 mg/mL (39.42 mM; Need ultrasonic) Storage Store at -80°C, protect from light

General For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the tips ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Evaluation sample solution : ship with blue ice All other available size: ship with RT, or Condition blue ice upon request.

Structure



Background

PROTAC CDK6 Degradar 1 (compound 48a) is a potent and selective PROTAC CDK6 degrader with a DC₅₀ of 0.037 μM. PROTAC CDK6 Degradar 1 exhibits selectivity over CDK4 (DC₅₀ > 10 μM). PROTAC CDK6 Degradar 1 induces G0/G1 cell-cycle arrest and apoptosis through inhibition of CDK6 downstream signaling. PROTAC CDK6 Degradar 1 reduces tumor burden and CDK6 levels in a MOLM-14 xenograft mouse model. PROTAC CDK6 Degradar 1 can be used for CDK6-driven cancers research, such as acute myeloid leukemia (AML) [1]. (Pink: CDK6 ligand; Blue: KLHDC2 ligand; Black: linker).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

[1]. Jeon E, et al. Selective CDK6 Degradation via the KLHDC2 E3 Ubiquitin Ligase. J Med Chem. 2025 Dec 11;68(23):25539-25568.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA