

Product Data Sheet

Product Name: SnapIt
Cat. No.: GA24161

Chemical Properties

Cas. No.

Formula M.Wt 197.23

Solubility Storage < -15°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Serine derivative containing a reactive moiety for chemoselective conjugation via inverse electron demand Diels-Alder reaction to tetrazines or 1,3-dipolar cycloaddition with azides. It reacts rapidly in water at physiological pH and room temperature. Norbornene also adds thiols (strain-promoted thiol-ene reaction).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA