

Product Data Sheet

Product Name: Spacer phosphoramidite 9

Cat. No.: GC81505

Chemical Properties

Cas. No. 146668-73-7


Formula $C_{36}H_{49}N_2O_7P$ M.Wt 652.76

Solubility Storage Store at -20°C, stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Spacer phosphoramidite 9 (Spacer-9 CEP) is a Phosphoramidite. Spacer phosphoramidite 9 can be used for oligonucleotide synthesis.

In Vitro, Spacer phosphoramidite 9 is a monomeric raw material that can be used for nucleic acid synthesis.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA