

Product Data Sheet

Product Name: (S,S) Bisoleoyl-lysobisphosphatidic acid

Cat. No.: GC81814

Chemical Properties

Cas. No. 288089-44-1

Formula $C_{42}H_{79}O_{10}P$ M.Wt 775.04

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

(S,S) Bisoleoyl-lysobisphosphatidic acid is a lysobisphosphatidic acid (LBPA). LBPA plays a role in modulating the cholesterol flux through endosomes.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA