

Product Data Sheet

Product Name: (S)-PHA533533

Cat. No.: GC27957

Chemical Properties

Cas. No. 437982-89-3

Chemical Name Benzeneacetamide, N-(5-cyclopropyl-1H-pyrazol-3-yl)- α -methyl-4-(2-oxo-1-pyrrolidinyl)-, (α S)-

SMILES C[C@H](C(NC1=NNC(C2CC2)=C1)=O)C3=CC=C(N4CCCC4=O)C=C3


Formula $C_{19}H_{22}N_4O_2$ M.Wt 338.40

Solubility DMSO : 100 mg/mL (295.51 mM; Need ultrasonic) Storage 4°C, protect from light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

(S)-PHA533533 is a CDK2 and CDK5 inhibitor with blood-brain barrier permeability, with IC₅₀ values of 37 nM and 55 nM, respectively. (S)-PHA533533 effectively restores UBE3A expression by downregulating UBE3A-ATS and relieving the epigenetic silencing of paternal UBE3A in mature neurons. (S)-PHA533533 can be used for the research of Angelman syndrome^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

In Vitro, (S)-PHA533533 (1 μ M; 72 h) potently activates paternal Ube3a expression in primary cortical neurons from Ube3a^{m+/pYFP} mice^[1].

(S)-PHA533533 (72 h) activates paternal Ube3a expression in primary cortical neurons from Ube3a^{m+/pYFP} mice, with an EC₅₀ of 0.54 μ M, an E_{MAX} of 80.02% paternal UBE3A-YFP-positive neurons, and a CC₅₀ of 2.98 μ M^[1].

(S)-PHA533533 (1 μ M; 72 h) activates paternal Ube3a expression in primary cortical neurons derived from Ube3a^{m+/pYFP} mice, and this process is independent of the inhibition of CDK2, CDK5 and TOP1^[1].

(S)-PHA533533 (72 h) downregulates UBE3A-ATS and relieves the silencing of paternal UBE3A in forebrain cortical neurons derived from ASdel1-0ΔI-P human induced pluripotent stem cells, with an EC₅₀ of 0.59 μ M for UBE3A induction and an IC₅₀ of 0.42 μ M for UBE3A-ATS downregulation^[1].

In Vivo, (S)-PHA533533 (2 mg/kg; i.p.; single injection) produces brain-wide unsilencing of paternal Ube3a, significantly increasing Ube3a mRNA and UBE3A protein levels in neurons of Angelman syndrome model mice^[1].

References:

[1]. Vihma H, et al. Ube3a unsilencer for the potential treatment of Angelman syndrome. Nat Commun. 2024;15(1):5558. Published 2024 Jul 8.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA