

Product Data Sheet

Product Name: JA310
Cat. No.: GC73860

Chemical Properties

Cas. No. 3022928-90-8


Formula $C_{17}H_{19}N_7O$ M.Wt 337.38

Solubility DMSO : 2 mg/mL (5.93 mM; ultrasonic and warming and heat to 60°C) Storage 4°C, stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

JA310 is a highly selective MST3 kinase inhibitor. JA310 has high cellular potency against MST3 with an EC50 value of 106 nM.

References:

[1]. Amrhein JA, et al. Synthesis of Pyrazole-Based Macrocycles Leads to a Highly Selective Inhibitor for MST3. J Med Chem. 2023 Dec 21.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA