

Product Data Sheet

Product Name: LY2922470

Cat. No.: GC36507

Chemical Properties

Cas. No. 1423018-12-5

SMILES COC1=C(N(CC2=CC=C(COC3=CC=C([C@@H](C#CC)CC(O)=O)C=C3)S2)CCC4)C4=CC=C1


Formula $C_{28}H_{29}NO_4S$ M.Wt 475.6

Solubility DMSO: ≥ 125 mg/mL (262.83 mM) Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility, please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

LY2922470 is a potent, selective and orally available agonist of the G protein-coupled receptor 40 (GPR40), with EC_{50} s of 7 nM, 1 nM and 3 nM for human GPR40, mouse GPR40 and rat GPR40, respectively. LY2922470 reduces glucose levels along with significant increases in insulin and GLP-1, is potential for the treatment of type 2 diabetes mellitus (T2DM)[1]. EC_{50} : 7 nM (human GPR40), 1 nM (mouse GPR40), 3 nM (rat GPR40)[1]

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. Hamdouchi C, et al. The Discovery, Preclinical, and Early Clinical Development of Potent and Selective GPR40 Agonists for the Treatment of Type 2 Diabetes Mellitus (LY2881835, LY2922083, and LY2922470). J Med Chem. 2016 Dec 22;59(24):10891-10916.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA