

Product Data Sheet

Product Name: Funapide (TV 45070)

Cat. No.: GC30628

Chemical Properties

Cas. No. 1259933-16-8

SMILES O=C([C@]12COC3=C1C=C(OCO4)C4=C3)N(CC5=CC=C(C(F)(F)F)O5)C6=C2C=CC=C6


Formula $C_{22}H_{14}F_3NO_5$ M.Wt 429.35

Solubility DMSO : 200 mg/mL (465.82 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Funapide (TV 45070; XEN402) is a potent Sodium Channel Nav1.7 inhibitor.

[1]. Price N , et al. Safety and Efficacy of a Topical Sodium Channel Inhibitor (TV-45070) in Patients With Postherpetic Neuralgia (PHN): A Randomized, Controlled, Proof-of-Concept, Crossover Study, With a Subgroup Analysis of the Nav1.7 R1150W Genotype. Clin J Pain. 2017 Apr;33(4):310-318.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA