

## Product Data Sheet

---

Product Name: Copper ionophore I

Cat. No.: GC81544

### Chemical Properties

Cas. No. 101264-90-8


Formula  $C_{13}H_{18}N_2O_2$  M.Wt 234.29

Solubility DMSO: 50 mg/mL (213.41 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



### Background

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

## Product Data Sheet

---

Copper ionophore I is an efficient copper ion carrier. Copper ionophore I can regulate various intracellular copper levels to induce cuproptosis, such as 4T1 (  $IC_{50} = 0.45 \mu M$ ) and MDA-MB-231 (  $IC_{50} = 1.21 \mu M$ ) cells. Copper ionophore I can induce an increase in ROS and cause mitochondrial dysfunction. Copper ionophore I can reduce the expression of FDX1, DLAT, and LIAS proteins. Copper ionophore I can activate immune related pathways and promote T cell infiltration. Copper ionophore I can be used for cancer research [1].

### References:

[1]. Yu L, et al. Rational design of copper ionophores for efficient induction of cuproptosis via simple n-alkyl modification. Eur J Med Chem. 2026 Jan 5;301:118257.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: [tech@glpbio.com](mailto:tech@glpbio.com)

Address: 10292 Central Ave. #205, Montclair, CA, USA