

## Product Data Sheet

Product Name: 5-(Piperidin-4-yl)-2-cyanophenylboronic acid

Cat. No.: GF23454

### Chemical Properties

Cas. No. 2225169-52-6

Formula C<sub>12</sub>H<sub>15</sub>BN<sub>2</sub>O<sub>2</sub>

M.Wt 230.0707

Solubility

Storage -20°C

**General tips** For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

**Shipping Condition** Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA