

Product Data Sheet

Product Name: CP 375
Cat. No.: GC30575

Chemical Properties

Cas. No. 752186-89-3

SMILES O=C1C(O)=C(C(OC)CC)N(C)C(C)=C1


Formula $C_{11}H_{17}NO_3$ M.Wt 211.26

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

CP 375 is a Fe³⁺ chelating agent, with a log K₁ value of 14.50.

CP 375 is a Fe³⁺ chelating agent, with a log K₁ value of 14.50, with the absolute deviation between the predicted and experimental values of 0.12 (14.62 vs. 14.50)[1].

[1]. Chen YL, et al. Prediction of 3-hydroxypyridin-4-one (HPO) log K₁ values for Fe(III). Dalton Trans. 2012 Sep 21;41(35):10784-91.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA