

Product Data Sheet

Product Name: Boc-Arg(Z)-OH

Cat. No.: GA11365

Chemical Properties

Cas. No. 51219-18-2

Chemical Name (2S)-5-[[amino(phenylmethoxycarbonylamino)methylidene]amino]-2-[(2-methylpropan-2-yl)oxycarbonylamino]pentanoic acid

SMILES CC(C)(C)OC(=O)NC(CCCN=C(N)NC(=O)OCC1=CC=CC=C1)C(=O)O

Formula $C_{19}H_{28}N_4O_6$

M.Wt 408.5


Solubility Soluble in water or 1% acetic acid

Storage Desiccate at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Boc-Arg(Z)-OH

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA