

Product Data Sheet

Product Name: PF-04979064

Cat. No.: GC19283

Chemical Properties

Cas. No. 1220699-06-8

SMILES O=C(N1C2CCN(C([C@@H](O)C)=O)CC2)N(C)C3=C1C4=NC(C5=CC=C(C)N=C5)=CC=C4N=C3


Formula $C_{24}H_{26}N_6O_3$ M.Wt 446.5

Solubility DMSO : 11 mg/mL (24.64 mM) Storage Store at -20°C

General For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the tips ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Evaluation sample solution : ship with blue ice All other available size: ship with RT, or Condition blue ice upon request.

Structure



Background

PF-04979064 is a potent and selective PI3K/mTOR dual kinase inhibitor with K_i of 0.13 nM and 1.42 nM for PI3K α and mTOR, respectively.

PF-04979064 is tested against human class I PI3K isoforms PI3K α , PI3K γ , and PI3K δ , with PI3K α K_i of 0.13 nM, PI3K γ K_i of 0.111 nM, and PI3K δ K_i of 0.122 nM[1].

PF-04979064 is progressed to rat in vivo PK studies and exhibits robust PK profile: V_{dss} =5.23 L/kg,

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Cl=19.3 mL/min/kg, T_{1/2}=1.85 h, and F%=61%[1].

References:

[1]. Cheng H et al. Discovery of the Highly Potent PI3K/mTOR Dual Inhibitor PF-04979064 through Structure-Based Drug Design. ACS Med Chem Lett, 2012 Nov 7, 4(1):91-7.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA