

Product Data Sheet

Product Name: CDK4/6-IN-11

Cat. No.: GC65190

Chemical Properties

Cas. No. 2139329-47-6


Formula C₄₃H₄₉N₁₁O₇ M.Wt 831.92

Solubility DMSO : 25 mg/mL (30.05 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

CDK4/6-IN-11 is a potent PROTAC CDK4/6 degrader.

[1]. Nathanael S. Gray, et al. Degradation of cyclin-dependent kinase 4/6 (cdk4/6) by conjugation of cdk4/6 inhibitors with e3 ligase ligand and methods of use. WO2020023480 A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA