

Product Data Sheet

Product Name: FSG67
Cat. No.: GC65176

Chemical Properties

Cas. No. 1158383-34-6


Formula C₁₆H₂₅NO₄S M.Wt 327.44

Solubility DMSO : 100 mg/mL (305.40 mM; Storage Store at -20°C
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

FSG67 is a glycerol 3-phosphate acyltransferase (GPAT) inhibitor with an IC₅₀ of 24 μM[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

FSG67 reduces GSK3 β phosphorylation and signaling after acetaminophen overdose[1]. FSG67 results in a dose-dependent increase in oxidative metabolism in mature adipocytes with an IC₅₀ of $27.7 \pm 4.4 \mu\text{M}$ [2].

[1]. Clemens MM, et al. The inhibitor of glycerol 3-phosphate acyltransferase FSG67 blunts liver regeneration after acetaminophen overdose by altering GSK3 β and Wnt/ β -catenin signaling. Food Chem Toxicol. 2019 Mar;125:279-288. [2]. Reilly SM, et al. Catecholamines suppress fatty acid re-esterification and increase oxidation in white adipocytes via STAT3. Nat Metab. 2020 Jul;2(7):620-634.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA