

Product Data Sheet

Product Name: Ac-Gly-Ala-Lys-AMC

Cat. No.: GA20567

Chemical Properties

Cas. No. 1926163-46-3

Formula $C_{23}H_{31}N_5O_6$ M.Wt 473.53

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Ac-GAK-AMC, control peptide for protease-coupled histone deacetylase (HDAC) assay with Ac-Gly-Ala-Lys(Ac)-AMC.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA