

Product Data Sheet

Product Name: DCAF1 ligand 2

Cat. No.: GC81017

Chemical Properties

Cas. No.


Formula $C_{31}H_{41}ClN_6O_2$ M.Wt 565.15

Solubility DMSO: 1 mg/mL (1.77 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

DCAF1 ligand 2 (Tert-butyl (2-((2-(1-(4-chlorophenyl)cyclohexyl)-7-(piperazin-1-yl)quinazolin-4-yl)amino)ethyl)carbamate) is a ligand of the E3 connecting enzyme DCAF1. DCAF1 ligand 2 can be used for the synthesis of PROTACs [1].

References:

[1]. Schröder M, et al. DCAF1-based PROTACs with activity against clinically validated targets overcoming intrinsic- and acquired-degrader resistance. Nat Commun. 2024 Jan 4;15(1):275. doi:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

10.1038/s41467-023-44237-4. Erratum in: Nat Commun. 2024 Aug 19;15(1):7109.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA