

Product Data Sheet

Product Name: Dilmapimod (SB-681323)

Cat. No.: GC31697

Chemical Properties

Cas. No. 444606-18-2

SMILES O=C1C=CC2=C(C3=CC=C(F)C=C3C)N=C(NC(CO)CO)N=C2N1C4=C(F)C=CC=C4F


Formula $C_{23}H_{19}F_3N_4O_3$ M.Wt 456.42

Solubility DMSO : 125 mg/mL (273.87 mM) Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Dilmapimod (SB-681323) is a potent p38 MAPK inhibitor that potentially suppresses inflammation in chronic obstructive pulmonary disease.

Dilmapimod (SB-681323) potently inhibits p38 MAPK signaling, measured by inhibition of sorbitol-induced pHSP27 expression[1].

[1]. Singh D, et al. A randomized, placebo-controlled study of the effects of the p38 MAPK inhibitor

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

SB-681323 on blood biomarkers of inflammation in COPD patients. J Clin Pharmacol. 2010 Jan;50(1):94-100.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA