

Product Data Sheet

Product Name: Timentin

Cat. No.: GF01034

Chemical Properties

Cas. No. 86482-18-0

Formula $C_{15}H_{16}N_2O_6S_2 \cdot C_8H_9NO_5$ M.Wt 583.593

Solubility DMSO : 100 mg/mL (171.35 mM; Storage -20°C
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Augpenin (Ticarcillin-clavulanic acid mixt.) is a mixture of Ticarcillin and Clavulanic acid in a ratio of 15:1.

[1]. Shiono A, et al. [Clinical evaluation of Augpenin (clavulanic acid-ticarcillin) in urinary tract infections: especially against beta-lactamase producing strains]. Hinyokika Kyo. 1993 Aug;39(8):777-82. Japanese.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA