

Product Data Sheet

Product Name: MAGE-3 Antigen (271-279) (human)

Cat. No.: GA23165

Chemical Properties

Cas. No. 160295-81-8

Formula $C_{53}H_{79}N_{13}O_{10}$

M.Wt 1058.29

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

This cytolytic T lymphocyte (CTL)-defined MAGE-3 antigen is associated with the human leukocyte antigen (HLA)-A2 molecule. The MAGE-3 gene is a member of a multigene family that is selectively expressed by subsets of different human tumor types, includ

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA