

Product Data Sheet

Product Name: TLQP-21 TFA

Cat. No.: GC61335

Chemical Properties

Cas. No.

Formula $C_{109}H_{171}F_3N_{40}O_{28}$ M.Wt 2546.77

Solubility H_2O : 50 mg/mL (19.63 mM; Need ultrasonic) Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

TLQPPASSRRRRHFHHALPPAR (TFA salt)

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Background

TLQP-21 TFA, a VGF-derived peptide endowed of endocrine and extraendocrine properties, is a potent G-protein-coupled receptor complement-3a receptor1 (C3aR1) agonist (EC50: mouse TLQP-21=10.3 μ M; human TLQP-21=68.8 μ M). TLQP-21 TFA activates C3aR1 to induce an increase of intracellular Ca²⁺. TLQP-21 TFA is used for the research in regulation of nociception and other relevant physiologic functions[1][2].

TLQP-21 TFA is a peptide of 21 amino acids. At a dose of 3 μ M TLQP-21 induces up to ~69% of the corresponding contraction promoted by acetylcholine[1][2].

[1]. Elena Bresciani , et al. TLQP-21, A VGF-Derived Peptide Endowed of Endocrine and Extraendocrine Properties: Focus on In Vitro Calcium Signaling. Int J Mol Sci. 2019 Dec 24;21(1):130. [2]. Cheryl Cero, et al. The TLQP-21 Peptide Activates the G-protein-coupled Receptor C3aR1 via a Folding-Upon-Binding Mechanism. ructure. 2014 Dec 2;22(12):1744-1753.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA