

Product Data Sheet

Product Name: (R)-MMG-1

Cat. No.: GC81484

Chemical Properties

Cas. No. 2064030-22-2

Formula $C_{35}H_{70}O_5$

M.Wt 570.93


Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

(R)-MMG-1 is the R-enantiomer of MMG-1. MMG-1, a monomycoloyl glycerol, is a C-type lectin receptor agonist. MMG-1, a lipid membrane stabilizer, is an immunopotentiator in liposomes. MMG-1 can be used to prepare lipid nanoparticles for drug delivery [1].

References:

[1]. Abhijeet G. Lokras, et al. Cholesterol in mRNA-Lipid Nanoparticles can be Replaced with the Synthetic Mycobacterial Monomycoloyl Glycerol Analogue MMG-1. Advanced Functional Materials. 15 August 2025.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA