

Product Data Sheet

Product Name: Cy3 labeled Inclisiran sodium

Cat. No.: GC80375

Chemical Properties

Cas. No.

Formula

M.Wt 16975.7

Solubility

Storage Store at -20°C, away from moisture

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Cy3 labeled Inclisiran sodium is a Cy3 labeled Inclisiran sodium.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA