

Product Data Sheet

Product Name: Diethylaminoethyl Sephadex A25

Cat. No.: GD05169

Chemical Properties

Cas. No. 12609-80-2

Formula C₉H₆FNO₂ M.Wt 0

Solubility DMF : < 1 mg/mL (insoluble); DMSO : < 1 mg/mL (insoluble or slightly soluble) Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Structure Not Available**

Background

DEAE Cross-linked dextran A 25 is a weak anion exchanger.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA