

Product Data Sheet

Product Name: AZD9496 deacrylic acid phenol

Cat. No.: GC74139

Chemical Properties

Cas. No. 2173404-70-9

Formula $C_{22}H_{23}F_3N_2O$


M.Wt 388.43

Solubility Storage -20°C, stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

AZD9496 deacrylic acid phenol is a fragment of target protein ligand of PROTAC ER Degradator-4 [1].

References:

[1]. Bin Yang, et al. Compounds and their use in treating cancer. WO2019123367A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA