

Product Data Sheet

Product Name: TB Dry
Cat. No.: GC13874

Chemical Properties

Cas. No.

Formula

M.Wt

Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

TB Dry is a single component high nutrient medium formulated to maximize yields of recombinant plasmids and cosmids. The most important advantage of TB DRY is that all components are combined into a single powder, eliminating the need to autoclave two separate solutions. TB Dry yields 4-7 times more plasmid DNA than LB broth.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA