

Product Data Sheet

Product Name: Fmoc-Phe(4-NO₂)-Wang resin

Cat. No.: GA10817

Chemical Properties

Cas. No. N/A

Formula C₂₄H₁₉N₂O₅Pol M.Wt 415.423

Solubility Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Fmoc-Phe(4-NO₂)-Wang resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA