

Product Data Sheet

Product Name: Ribosomal protein L3 peptide (202-222) amide

Cat. No.: GP10147

Chemical Properties

Cas. No.


Formula $C_{114}H_{182}N_{42}O_{25}S$ M.Wt 2573

Solubility $\geq 257.3\text{mg/mL}$ in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA