

Product Data Sheet

Product Name: 3-(Piperidin-4-yl)-5-(iso-butyl)phenylboronic acid

Cat. No.: GF32993

Chemical Properties

Cas. No. 2225174-85-4

Formula C₁₅H₂₄BN₂O₂

M.Wt 261.16756

Solubility

Storage -20°C, sealed storage, away from moisture

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA