

Product Data Sheet

Product Name: L2H2-6OTD intermediate-2

Cat. No.: GC74148

Chemical Properties

Cas. No. 890403-75-5

Formula $C_{14}H_{20}N_2O_6$

M.Wt 312.32


Solubility

Storage 4°C, stored under nitrogen, away from moisture

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

L2H2-6OTD intermediate-2 is an intermediate of telomerase inhibitor, which can be used in the synthesis of ADC.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

[1]. Kazuaki Shibata, et al. Synthesis of Heptaoxazole Macrocyclic Analogues of Telomestatin and Evaluation of Their Telomerase Inhibitory Activities. Bulletin of the Chemical Society of Japan. Volume 86, Issue 12, December 2013, Pages 1453-1465.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA