

Product Data Sheet

Product Name: N-Boc-4-piperidinemethanol

Cat. No.: GD02107

Chemical Properties

Cas. No. 123855-51-6

Formula C₁₁H₂₁NO₃

M.Wt 215.289


Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

N-Boc-4-piperidinemethanol (compound 4) is an acetylcholinesterase (AChE) inhibitor. N-Boc-4-piperidinemethanol can be used to study Alzheimer's Disease.

[1]. Villalobos A, et al. Novel benzisoxazole derivatives as potent and selective inhibitors of acetylcholinesterase[J]. Journal of medicinal chemistry, 1994, 37(17): 2721-2734.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA