

Product Data Sheet

Product Name: Fmoc-Lys(Tfa)-OH

Cat. No.: GA10918

Chemical Properties

Cas. No. 76265-69-5

Chemical Name (2S)-2-(9H-fluoren-9-ylmethoxycarbonylamino)-6-[(2,2,2-trifluoroacetyl)amino]hexanoic acid

SMILES C1=CC=C2C(=C1)C(C3=CC=CC=C32)COC(=O)NC(CCCCNC(=O)C(F)(F)F)C(=O)O

Formula $C_{23}H_{23}F_3N_2O_5$

M.Wt 464.4


Solubility

Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Fmoc-Lys(Tfa)-OH

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA