

Product Data Sheet

Product Name: HO-PEG24-OH

Cat. No.: GC63915

Chemical Properties

Cas. No. 2243942-52-9

Formula C₄₈H₉₈O₂₅ M.Wt 1075.28

Solubility DMSO : 50 mg/mL (46.50 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Structure



www.glpbio.com