

Product Data Sheet

Product Name: Thalidomide-piperazine-Boc

Cat. No.: GC63220

Chemical Properties

Cas. No. 2222114-64-7


Formula $C_{22}H_{26}N_4O_6$ M.Wt 442.47

Solubility DMSO : 41.67 mg/mL (94.18 mM; Storage Store at -20°C
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Thalidomide-piperazine-Boc is an intermediate that can be used in the synthesis of B-cell lymphoma 6 protein (BCL6) PROTAC[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. Michael Berlin, et al. Bifunctional molecules containing an e3 ubiquitine ligase binding moiety linked to a bcl6 targeting moiety. WO2021077010A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA