

## Product Data Sheet

Product Name: Isoxazole

Cat. No.: GD15157

### Chemical Properties

Cas. No. 288-14-2

Formula

M.Wt

Solubility DMSO : 100 mg/mL (1448.02 mM; Need ultrasonic)


Storage

2-8°C argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



### Background

Isoxazole is a agent scaffold that can be used to develop new agents with various biological activities (e.g. anti-microbial, anti-viral, anti-cancer, anti-inflammatory, immunomodulatory, anti-convulsant or anti-diabetic properties).

[1]. Zhu J, et al. The recent progress of isoxazole in medicinal chemistry. Bioorg Med Chem. 2018 Jul 23;26(12):3065-3075.

[2]. Sysak A, et al. Isoxazole ring as a useful scaffold in a search for new therapeutic agents. Eur J

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

## Product Data Sheet

---

Med Chem. 2017 Sep 8;137:292-309.

***Caution: Product has not been fully validated for medical applications. For research use only.***

***Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: [tech@glpbio.com](mailto:tech@glpbio.com)***

***Address: 10292 Central Ave. #205, Montclair, CA, USA***