

Product Data Sheet

Product Name: ARM1
Cat. No.: GC11201

Chemical Properties

Cas. No. 68729-05-5

Chemical Name 4-[4-(phenylmethyl)phenyl]-2-thiazolamine

SMILES NC(N1)=SC=C1C(C=C2)=CC=C2CC3=CC=CC=C3


Formula $C_{16}H_{16}N_2S$ M.Wt 268.4

Solubility ≤ 1 mg/ml in ethanol; 10mg/ml in DMSO; 15mg/ml in dimethyl formamide
Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

ARM1 (4-(4-benzylphenyl)thiazol-2-amine) is a LTB₄ synthesis inhibitor [1].

Leukotriene A₄ (LTA₄) hydrolase/aminopeptidase is a bifunctional zinc metalloenzyme that catalyzes the synthesis of the proinflammatory mediator LTB₄ from LTA₄ and cleaves the chemotactic

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

tripeptide Pro-Gly-Pro. LTB₄ is a key lipid mediator in the innate immune response [1][2].

ARM1 is a LTB₄ synthesis inhibitor. ARM1 is a thiazolamine that inhibits LTB₄ synthesis and conversion of LTA₄ to LTB₄ by purified LTA₄H with K_i value of 2.3 μ M. 50- to 100-fold higher concentrations of ARM1 did not significantly affect hydrolysis of Pro-Gly-Pro by LTA₄ hydrolase. In isolated human polymorphonuclear neutrophils (PMNs), ARM1 efficiently inhibited LTB₄ synthesis with IC₅₀ value of about 0.5 μ M and complete inhibition at 5 μ M. ARM1 inhibited the epoxide hydrolase activity of LTA₄H without significantly affecting its ability to cleave Pro-Gly-Pro. In the crystal structure of LTA₄H in complex with ARM1, ARM1 occupied the binding site for the ω -end of LTA₄ [1].

References:

- [1]. Stsiapanava A, Olsson U, Wan M, et al. Binding of Pro-Gly-Pro at the active site of leukotriene A₄ hydrolase/aminopeptidase and development of an epoxide hydrolase selective inhibitor. Proc Natl Acad Sci U S A. 2014 Mar 18;111(11):4227-32.
- [2]. Haeggström JZ. Leukotriene A₄ hydrolase/aminopeptidase, the gatekeeper of chemotactic leukotriene B₄ biosynthesis. J Biol Chem. 2004 Dec 3;279(49):50639-42.
- [3]. Snelgrove RJ, Jackson PL, Hardison MT, et al. A critical role for LTA₄H in limiting chronic pulmonary neutrophilic inflammation. Science. 2010 Oct 1;330(6000):90-4.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA