

Product Data Sheet

Product Name: HCy-Lyso

Cat. No.: GC91751

Chemical Properties

Cas. No. N/A


Formula $C_{24}H_{30}N_2O$ M.Wt 362.5

Solubility Chloroform: Sparingly Soluble: 1-10 mg/ml Storage $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

HCy-Lyso is a lysosome-specific fluorescent reporter for the detection of hydroxyl radicals. [1] It selectively fluoresces in low pH environments in the presence of hydroxy radicals over hypochlorite, superoxide, nitric oxide, or peroxynitrite radicals. HCy-Lyso ($10\text{ }\mu\text{M}$) displays lysosomal fluorescence in the presence of the ferroptosis inducers erastin or (1S,3R)-RSL3 , an effect that can be blocked by the ferroptosis inhibitor ferrostatin-1 , in 4T1 cells. HCy-Lyso displays excitation/emission maxima of 510/592 nm, respectively, in low pH environments in the presence of hydroxy radicals.

References:

[1].Zhong, L., Fu, D., Xu, J., et al. Rational design of a lysosome-targeted fluorescent probe for monitoring the generation of hydroxyl radicals in ferroptosis pathways RSC Adv. 14(18)12864-12872(2024).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA