

Product Name: (S,R,S)-AHPC-Me-NHBoc
Cat. No.: GC81015

Cas. No. 1997302-16-5


Formula	C ₂₈ H ₄₀ N ₄ O ₅ S	M.Wt	544.71
---------	---	------	--------

Solubility	DMSO: 50 mg/mL (91.79 mM; Need ultrasonic)	Storage	Store at -20°C, stored under nitrogen
------------	--	---------	---------------------------------------

General tips	For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.
--------------	---

Shipping Condition	Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.
--------------------	--

Structure



(S,R,S)-AHPC-Me-NHBoc is an E3 ubiquitin ligase VHL ligand used to recruit the VHL protein.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

(S,R,S)-AHPC-Me-NHBoc can be linked to a target protein ligand via a linker to form a PROTAC.

In Vitro, PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA