

Product Data Sheet

Product Name: Tecadenoson (CVT-510)

Cat. No.: GC30851

Chemical Properties

Cas. No. 204512-90-3

SMILES OC[C@@H]1[C@H]([C@H]([C@H](N2C=NC3=C2N=CN=C3N[C@H]4COCC4)O1)O)O


Formula $C_{14}H_{19}N_5O_5$ M.Wt 337.33

Solubility DMSO : ≥ 155 mg/mL (459.49 mM) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Protocol

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

**Kinase
experiment:**

The effect of Tecadenoson on binding to A1 and A2A-adenosine receptors of porcine forebrain and striatum membranes, respectively, are determined. Assays for A1 and A2A receptors are carried out by using the A1 receptor antagonist [3H]CPX and the A2A receptor agonist [3H]CGS 21680. Membranes are treated with adenosine deaminase (2 U/mL) for 20 min at room temperature prior to and during radioligand binding assays. Membranes (0.2-0.7 mg), adenosine deaminase, and the indicated radioligand are incubated for 3 h in a 300 μ L volume of Tris-HCl buffer (50 mM) (pH 7.4). Assays are carried out in triplicate at room temperature. After the incubation period, bound and free radioligand are diluted by the addition of ice-cold Tris-HCl buffer (5 mL), and immediately separated by vacuum filtration of assay contents onto Whatman GF/C filters and washing of trapped membranes with Tris-HCl buffer (20 mL). Filter disks containing membrane-bound radioactivity are placed in 4 mL Scintiverse, and the radioactivity is quantified by a liquid scintillation counter. Specific binding of [3H]CPX and [3H]CGS 21680 is defined as membrane binding displaced in the presence of CPT (10 μ M) and R-PIA (10 μ M), respectively[1].

**Animal
experiment:**

Rat: The effects of Tecadenoson on heart rate and to reduce serum NEFA concentration are determined in separate groups of rats to avoid the effects of animal handling and blood sampling on heart rate. Three days before an experiment, a catheter (0.025-mm outer diameter) is implanted in the left common carotid artery of each rat using aseptic conditions and sterile technique. The catheter is tunneled subcutaneously to the dorsal surface. After recovery from anesthesia, rats are placed in metabolic cages to facilitate handling and blood sampling. Blood samples (0.2 mL) are drawn before and at various time points after i.p. injection of either Tecadenoson or vehicle (DMSO in saline). A 0.4-mL volume of 1% sodium citrate in saline is administered after withdrawal of each blood sample to replace blood volume and prevent clotting in the carotid artery catheter. Serum is collected from each sample after centrifugation of the clotted blood. Serum samples are stored at -80°C until analysis. Serum NEFA concentration is determined using an enzymatic colorimetric assay kit[2].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

- [1]. Snowdy S, et al. A comparison of an A1 adenosine receptor agonist (CVT-510) with diltiazem for slowing of AVnodal conduction in guinea-pig. Br J Pharmacol. 1999 Jan;126(1):137-46.
- [2]. Fraser H, et al. N-[3-(R)-tetrahydrofuranyl]-6-aminopurine riboside, an A1 adenosine receptor agonist, antagonizes catecholamine-induced lipolysis without cardiovascular effects in awake rats. J Pharmacol Exp Ther. 2003 Apr;305(1):225-31.

Background

Tecadenoson (CVT-510) is a selective A1 adenosine receptor agonist.

In the atrial-paced isolated heart, Tecadenoson is approximately 5 fold more potent to prolong the stimulus-to-His bundle (S-H interval), a measure of slowing AV nodal conduction ($EC_{50}=41$ nM) than to increase coronary conductance ($EC_{50}=200$ nM). At concentrations of Tecadenoson (40 nM) and diltiazem (1 μ M) that causes equal prolongation of S-H interval (~ 10 ms), diltiazem, but not Tecadenoson, significantly reduces left ventricular developed pressure (LVP) and markedly increases coronary conductance. Tecadenoson shortens atrial ($EC_{50}=73$ nM) but not the ventricular monophasic action potentials (MAP)[1].

In atrial-paced anaesthetized guinea-pigs, intravenous infusions of Tecadenoson and diltiazem causes nearly equal prolongations of P-R interval[1]. Tecadenoson (2, 5, 20 μ g/kg i.p.) causes a rapid and sustained dose-dependent decrease in NEFA at doses that do not cause bradycardia. Tecadenoson given at 50 μ g/kg causes a significant bradycardia (50% decrease in heart rate at 25 min[2].

- [1]. Snowdy S, et al. A comparison of an A1 adenosine receptor agonist (CVT-510) with diltiazem for slowing of AVnodal conduction in guinea-pig. Br J Pharmacol. 1999 Jan;126(1):137-46. [2]. Fraser H, et al. N-[3-(R)-tetrahydrofuranyl]-6-aminopurine riboside, an A1 adenosine receptor agonist, antagonizes catecholamine-induced lipolysis without cardiovascular effects in awake rats. J Pharmacol Exp Ther. 2003 Apr;305(1):225-31.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA