

Product Data Sheet

Product Name: Zeclantamig

Cat. No.: GC79350

Chemical Properties

Cas. No. 3052965-90-6

Formula M.Wt 125.88

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Zeclantamig

Background

Zeclantamig is an immunoglobuline (H-gamma1_L-kappa)_scFv_{hk}-G1(h-CH₂-CH₃), anti-[Homo sapiens EPCAM et anti-[Homo sapiens CD3E], anticorps monoclonal chimérique, bispécifique, trivalent [1].

References:

[1]. WHO Drug Information, Vol. 39, No. 4, 2025

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA