

Product Data Sheet

Product Name: CP05
Cat. No.: GC81554

Chemical Properties

Cas. No. 1824650-30-7


Formula $C_{56}H_{99}N_{21}O_{18}S_2$ M.Wt 1418.65

Solubility DMSO: 100 mg/mL (70.49 mM; Need Storage Sealed storage, away from moisture ultrasonic)

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

CP05 is a peptide that specifically recognizes an exosome surface marker (CD63). CP05 anchors exosomes to target tissues for in vivo delivery of exosomes [1].

References:

[1]. Li X, et al. A drug delivery system constructed by a fusion peptide capturing exosomes targets to

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

titanium implants accurately resulting the enhancement of osseointegration peri-implant. Biomater Res. 2022 Dec 27;26(1):89.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA