

Product Data Sheet

Product Name: CleanAmp dUTP Set

Cat. No.: GN20177

Chemical Properties

Purity	>98%	Extinction Coefficient	
Formula		M.Wt	
Salt Form		Concentration	
Buffer		Storage	-20°C or below
Synonyms		Backbone	Phosphodiester
Base Analog	Adenosine, Cytidine, Guanosine, Uridine	Sugar Type	
Nucleotide Category			
Applications	In vitro Diagnostic (IVD), Laboratory Developed Tests (LDT) & Analyte Specific Reagents (ASR), PCR & RT-PCR		

Background

1 (50 mM) vial each of dATP, dCTP, dGTP and dUTP Each dNTP is analyzed separately by AX-HPLC, RP-HPLC, 31P NMR and Mass Spec. The components are tested in a standardized PCR assay.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA