

Product Data Sheet

Product Name: DYRK2-IN-2

Cat. No.: GC81307

Chemical Properties

Cas. No. 2416195-65-6

Formula $C_{21}H_{24}N_2O_3S$

M.Wt 384.49


Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

DYRK2-IN-2 (Compound C17) is a selective inhibitor of DYRK2, with its IC_{50} value being 40.3 nM. The inhibitory activity of DYRK2-IN-2 on DYRK1A (IC_{50} = 1842 nM), DYRK1B (IC_{50} = 1335 nM), DYRK4 (IC_{50} = 1931 nM), DYRK3 (IC_{50} = 112 nM), and CLKs (IC_{50} = 540-6496 nM) is relatively low. DYRK2-IN-2 inhibits the phosphorylation of Tau protein at Thr212 and shows moderate cytotoxicity in HT22 cells. DYRK2-IN-2 can be used in cancer research [1].

References:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. Lindberg MF, et al. Comparative Efficacy and Selectivity of Pharmacological Inhibitors of DYRK and CLK Protein Kinases. J Med Chem. 2023 Mar 23;66(6):4106-4130.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA