

Product Data Sheet

Product Name: GMC 2-29

Cat. No.: GB40218

Chemical Properties

Cas. No. 148672-15-5

Formula $C_{29}H_{31}N_5O_3$

M.Wt 497.59

Solubility

Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Selective 5-HT1B/1D antagonist; GR 127935-like drug with even greater potency in function assays

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA