

Product Data Sheet

Product Name: GPR55 agonist 4

Cat. No.: GC73616

Chemical Properties

Cas. No. 3006105-55-8


Formula $C_{19}H_{16}FN_5O_2$ M.Wt 365.36

Solubility DMSO : 100 mg/mL (273.70 mM; ultrasonic and warming and heat to 60°C) Storage -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

GPR55 agonist 4 (Compound 28) is a GPR55 agonist (EC₅₀: 131 nM, and 1.41 nM for hGPR55 and rGPR55). GPR55 agonist 4 induces β -arrestin recruitment to human GPR55.

References:

[1]. Hewer RC, et al. Discovery and Characterization of Novel CNS-Penetrant GPR55 Agonists. J Med Chem. 2023 Sep 28;66(18):12858-12876.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA