

Product Data Sheet

Product Name: H-Ile-2-Chlorotrityl Resin

Cat. No.: GA10732

Chemical Properties

Cas. No. N/A

Chemical Name (2S,3S)-2-amino-3-methylpentanoic acid

SMILES O=C([C@@]([C@](C([H])([H])[H])([H])C([H])([H])C([H])([H])[H])([H])N([H])[H])O[H]


Formula 100-200:mesh,1%DVB M.Wt 0.3-0.8mmol/g

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

H-Ile-2-Chlorotrityl Resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA