

Product Data Sheet

Product Name: DL-m-Tyrosine

Cat. No.: GA10798

Chemical Properties

Cas. No. 775-06-4

Chemical Name 2-amino-3-(3-hydroxyphenyl)propanoic acid

SMILES C1=CC(=CC(=C1)O)CC(C(=O)O)N


Formula $C_9H_{11}NO_3$ M.Wt 181.2

Solubility H_2O : 1.43 mg/mL (7.89 mM; Need ultrasonic) Storage Store at 2-8°C, protect from light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

DL-m-Tyrosine

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA