

Product Data Sheet

Product Name: ZNL-05-044

Cat. No.: GC73597

Chemical Properties

Cas. No.


Formula $C_{21}H_{22}Cl_2N_6OS$ M.Wt 477.41

Solubility DMSO : 100 mg/mL (209.46 mM; Storage 4°C, stored under nitrogen
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

ZNL-05-044 is a CDK11 inhibitor with an IC₅₀s of 0.23 μ M and 0.27 μ M against CDK11A and CDK11B, respectively (NanoBRET assay). ZNL-05-044 leads to G2/M cell cycle arrest and impairs RNA splicing.

References:

[1]. Li Z, et al. Synthesis and Structure-Activity relationships of cyclin-dependent kinase 11 inhibitors based on a diaminothiazole scaffold. Eur J Med Chem. 2022 Aug 5;238:114433.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA