

Product Data Sheet

Product Name: cis-KIN-8194

Cat. No.: GC73777

Chemical Properties

Cas. No. 330792-60-4

Formula C₂₈H₃₃N₇O

M.Wt 483.61


Solubility

Storage -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

cis-KIN-8194 is the cis-isomer of KIN-8194

References:

[1]. Lantermans HC, et al. The dual HCK/BTK inhibitor KIN-8194 impairs growth and integrin-mediated adhesion of BTKi-resistant mantle cell lymphoma. Leukemia. 2024 Mar 7.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA