

Product Data Sheet

Product Name: Fmoc-Cys(Acm)-Wang resin

Cat. No.: GA11063

Chemical Properties

Cas. No. N/A

Chemical Name (Z)-N-((((R)-2-((E)-(((9H-fluoren-9-yl)methoxy)(hydroxy)methylene)amino)-3-oxopropyl)thio)methyl)acetimidic acid

SMILES C/C(O)=N/CSC[C@@H](/N=C(OCC1C2=CC=CC=C2C3=CC=CC=C31)\O)([H])C=O


Formula $C_{21}H_{22}N_2O_4S$ M.Wt 398.48

Solubility Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Fmoc-Cys(Acm)-Wang resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA